



AREA of DETAIL



Adduct formation

Site of Aggregation



N-terminal pyroE formation

Conjugation Site (ADC)

Conjugation Site (ADC)

ThermoFisher
S C I E N T I F I C

Future analytical trends for the analysis of biotherapeutic proteins

Immunoglobulin protein | ca. 150,000 Daltons | participates in the immune response by acting as the antibody for a specific antigen | There are five main types: IgA, IgG, IgM, IgE, and IgD

Simon Cubbon
Pharma and Biopharma Team
Thermo Fisher Scientific


Humanized IgG antibody fragment (Fab) | 50,000 Daltons | VH, CH1 and VL, CL regions, linked by an intramolecular disulfide bond.

STRUCTURAL INSIGHTS

Industry Challenges

- Regulatory landscape is changing
- Speed, flexibility, quality and cost control: getting it right the first time and every time
- Driving towards continuous production – more production analysis
- Simple robust workflows required
- Proving bioequivalence and monitoring batch-to-batch variation
- Increasing biotherapeutic complexity is forcing new methodologies

A Complex Problem: Drug Safety and Quality



Safety

Is the product safe to use?
(e.g. Immunogenic effects?)

Potency

Does the drug have the expected effect?
(e.g. CDR complementation)

Knowledge

How do changes effect the therapeutic?
(e.g. Oxidation)

Quality

How do changes in process effect the product?
(e.g. Glucose concentration on glycoforms)


A Complex Problem: Drug Safety and Quality

- Can we reduce the analytical complexity?
- Can we get this data from fewer methods?
- Can these new methods be simple and robust enough for QC or even Process monitoring?

Multi Attribute Methods

One injection?


Sub-unit Analysis




Glycan Analysis



Released Glycan Analysis




Peptide Mapping




Intact Analysis



Aggregate Screening




Charge Variant Screening



Peptide Mapping: A Complete Multi-Attribute Method



Method Evaluation and
Advancements



Digestion

SMART Digest Kits



Separation

Acclaim Vanquish
Vanquish UHPLC




Detection

Exactive or Q Exactive Series



Processing

Chromeleon
BioPharma Finder



SMART

- Really easy to use
- Highly reproducible
- Rapid digestions



Separation

- Class leading, biocompatible
- Longer column lifetimes
- Unrivalled reproducibility



Detection




- Walk-up simplicity
- Minimal instrument contact
- Rock-steady results



Processing

- Easy to use and learn
- Enterprise & LIMS ready
- Compatible with other vendors instrumentation

Thermo Scientific SMART Digest Kits and Automation




Denaturation

Digestion

Clean up

5 different Rituximab Digests by 5 different Seminar Attendees

4.5 6.0 7.0 8.0 9.0 10.0 11.0 12.0 13.0 14.0 15.0 16.0 16.9




| Peak 1 | Peak 2 | Peak 3 | Peak 4 | Peak 5 | Peak 6 | Peak 7 | Peak 8 | Peak 9 | Peak 10 |
|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| 2.54 | 2.41 | 1.89 | 3.39 | 3.53 | 2.16 | 4.41 | 2.10 | 2.10 | 3.65 |
| Peak 11 | Peak 12 | Peak 13 | Peak 14 | Peak 15 | Peak 16 | Peak 17 | Peak 18 | Peak 19 | Peak 20 |
| 1.96 | 3.51 | 3.72 | 2.26 | 2.91 | 1.97 | 3.28 | 2.62 | 3.16 | 1.20 |

Advantages of MAM

| Required CQA Characterizations | Covered by SEC MAM? | CEX | rCE-SDS | nrCE-SDS | HILIC | ID ELISA | HCP ELISA |
|--------------------------------|---------------------|----------|---------|----------|-------|----------|-----------|
| CDR Tryptophan Degradation | Yes | No | No | No | No | No | No |
| C-terminal Amidation | Yes | Indirect | No | No | No | No | No |
| C-terminal Lysine | Yes | Yes | No | No | No | No | No |
| Cysteine Adducts | Maybe | Maybe | No | No | No | No | No |
| Deamidation | Yes | Indirect | No | No | No | No | No |
| Disulfide Isoforms | Maybe | Indirect | No | Yes | No | No | No |
| Disulfide Reduction | Maybe | No | No | Yes | No | No | No |
| Fragmentation (Peptide Bond) | Maybe | No | Yes | Yes | No | No | No |
| Fucosylation | Yes | No | No | No | No | No | No |
| Galactosylation | Yes | No | No | No | No | No | No |
| Glycation | Yes | No | Yes | Yes | No | No | No |
| HCP | Yes | No | No | No | No | No | Yes |
| High Mannose | Yes | No | No | No | Yes | No | No |
| Hydroxylysine | Yes | No | No | No | No | No | No |
| Identity | Yes | Yes | No | No | No | Yes | No |
| Methionine Oxidation | Yes | No | No | No | No | No | No |
| Mutations & Misincorporations | Yes | No | No | No | No | No | No |
| Non-concensus Glycosylation | Yes | No | Maybe | Maybe | No | No | No |
| Non-glycosylated Heavy Chain | Yes | No | No | No | No | No | No |
| N-terminal pyroGlutamate | Yes | Indirect | No | No | No | No | No |
| O-linked Glycans | Maybe | No | No | No | No | No | No |
| Residual Protein A | Yes | No | No | No | No | No | No |
| Signal Peptide | Yes | No | No | No | No | No | No |
| Thioether | Yes | No | No | No | No | No | No |
| Trisulfide | Maybe | No | No | No | No | No | No |
| Unusual Glycosylation | Yes | Indirect | Maybe | Maybe | Yes | No | No |

Building Targeted List of Critical Quality Attributes




CQA Overview

Sequence Details

| | | | |
|----------------------------|-------------------------|-------------|--------------------|
| Name: | HMSVictory | Created On: | 09/Nov/16 14:49:34 |
| Instrument: | QEHF_TheGhost | Updated On: | 20/Dec/16 17:07:26 |
| Imported Data: | False | | |
| First Injection: | NISTmAb_Control1_DD...1 | | |
| Processing Method: | SuperMAM | | |
| MS Acquisition Time [min]: | 70.51 | | |
| Method Length [min]: | 115.00 | | |
| Total Time [hrs]: | 63.25 | | |
| Data Vault: | ChromleonLocal | Created By: | Thermo |
| No. of Injections: | 33 | Updated By: | Thermo |


Deamidations and Isomerizations

| No. | Name | PENNKY Deam | FNWY IsoD | FNWY Deam | NOVWLK Deam | VVSV Deam | WQQ Deam | YCAR Deam | VDNAL Deam | SGTAS Deam |
|-----|-----------------------|-------------|-----------|-----------|-------------|-----------|----------|-----------|------------|------------|
| | | HC N391 | HC D283 | HC N283 | HC N78 | HC N317 | HC N436 | HC N85 | LC N151 | LC N136 |
| 11 | Control_Digest1_Rep1 | 1.58 | 0.11 | 0.04 | 0.03 | 0.05 | 0.23 | 0.00 | 0.21 | 8.3 |
| 12 | Control_Digest1_Rep2 | 1.59 | 0.11 | 0.04 | 0.03 | 0.05 | 0.24 | 0.00 | 0.30 | 7.9 |
| 13 | Control_Digest1_Rep3 | 1.58 | 0.11 | 0.03 | 0.03 | 0.05 | 0.24 | 0.00 | 0.32 | 8.5 |
| 15 | Control_Digest2_Rep1 | 1.54 | 0.11 | 0.04 | 0.04 | 0.05 | 0.23 | 0.00 | 0.40 | 7.9 |
| 18 | Control_Digest3_Rep1 | 1.48 | 0.10 | 0.04 | 0.04 | 0.04 | 0.40 | 0.00 | 0.38 | 8.4 |
| 24 | Stressed_Digest2_Rep1 | 10.89 | 22.34 | 6.66 | 0.04 | 11.78 | 6.80 | 0.12 | 1.94 | 31.6 |
| 25 | Stressed_Digest2_Rep2 | 11.03 | 22.14 | 6.55 | 0.04 | 11.98 | 6.80 | 0.12 | 1.68 | 31.1 |
| 26 | Stressed_Digest2_Rep3 | 11.25 | 22.41 | 6.31 | 0.04 | 11.45 | 6.84 | 0.12 | 1.65 | 30.7 |




Oxidations

| No. | Injection Name | DTLM Oxid | DIQM Oxid | DMIF Oxid | FNWY Oxid |
|-----|-----------------------|-----------|-----------|-----------|-----------|
| | | HC M253 | LC M4 | HC M100 | HC W279 |
| 11 | Control_Digest1_Rep1 | 1.4 | 0.3 | 1.0 | 0.01 |
| 12 | Control_Digest1_Rep2 | 1.4 | 0.3 | 1.0 | 0.01 |
| 13 | Control_Digest1_Rep3 | 1.5 | 0.3 | 1.1 | 0.01 |
| 15 | Control_Digest2_Rep1 | 1.3 | 0.2 | 0.7 | 0.01 |
| 18 | Control_Digest3_Rep1 | 1.2 | 0.2 | 0.7 | 0.01 |
| 24 | Stressed_Digest2_Rep1 | 2.9 | 0.2 | 1.1 | 0.02 |
| 25 | Stressed_Digest2_Rep2 | 2.9 | 0.2 | 1.2 | 0.02 |
| 26 | Stressed_Digest2_Rep3 | 3.0 | 0.2 | 1.2 | 0.02 |



Lysine Modifications



| No. | Injection Name | SLSSPGK | VSNK | APK |
|-----|-----------------------|---------------|-------------|------------|
| | | Lysine Trunc. | HC K329 Gly | LC K44 Gly |
| 11 | Control_Digest1_Rep1 | 28.4 | 0.6 | 0.1 |
| 12 | Control_Digest1_Rep2 | 28.2 | 0.5 | 0.1 |
| 13 | Control_Digest1_Rep3 | 27.2 | 0.6 | 0.1 |
| 18 | Control_Digest3_Rep1 | 26.6 | 0.6 | 0.2 |
| 24 | Stressed_Digest2_Rep1 | 25.6 | 0.3 | 0.1 |
| 25 | Stressed_Digest2_Rep2 | 25.9 | 0.3 | 0.1 |
| 26 | Stressed_Digest2_Rep3 | 25.5 | 0.3 | 0.1 |



Custom MS
Report Template

pH gradient elution


- Based on pI of protein
- Loss of retention with progressing pH gradient, depending on pI
- “Single” binding event, trapping at $pH < pI$ (for CEX)




Isoelectric Focusing on a Cation Exchange Column

- mAb binds to cation exchange sites on the column
- A gradient of increasing pH is applied
- mAb is released from the exchange site when the net charge on the mAb is neutral
- This interaction happens once, then the mAb runs through the rest of the column
- Column length has little effect on the resolution
- This is a concentration technique

Trastuzumab pH Gradient with Proprietary Volatile eluents



Direct MS analysis of Trastuzumab by Ion Exchange / Native Intact MS




| Glycoform | Theoretical Av. Mass | Experimental Av. Mass | D Mass (ppm) |
|-------------------|----------------------|-----------------------|--------------|
| (G0F)2 | 148056.2 | 148055.7 | 3.4 |
| G0F/G1F | 148218.3 | 148217.6 | 5.0 |
| G0F/G2F or (G1F)2 | 148380.5 | 148380.7 | -1.5 |
| G1F/G2F | 148542.6 | 148543.2 | -3.9 |




Conclusions


- Automation of protein digestion for precision and reproducibility in a robust peptide mapping workflow to introduce a MAM
- Volatile pH gradient elution of proteins for CVA has several advantages; global applicability, high loading capacity, easy method transfer, native MS compatibility, simple sample preparation
- Simple parent ion MS at high resolution, allowing on-line positive identifications at intact and peptide level
- One Injection, reliable and robust: several critical quality attributes
- New workflows for the future enabling characterization of several attributes in one injection provides ease of use and time saving



Peptide Mapping



Charge Variant Screening



Aggregate Screening



Intact Analysis



Glycan Analysis