The Power of Single Molecule Counting: The Future of Immunoassays

Sarah J. Hamren¹, Robert J. Goodson¹, Dr. Leslie Gordon² ¹Millipore Sigma, 25801 Industrial Blvd, Hayward, CA 94545, ²Progeria Research Foundation, 200 Lake St. Peabody, MA 01960

Introduction

The increasing need for new therapeutics that target different mechanisms of action, and methods for earlier diagnosis require increased understanding of disease processes and progression. Fundamental to this is the ability to distinguish between healthy and disease states. The clinical use of protein biomarkers often necessitates the measurement of small changes at low concentrations, in small volumes, which demands technology that can provide sensitive, precise and accurate quantification. Since the introduction of immunoassay (IA) technology and its application to the measurement of proteins in biological samples, significant improvements have been made in performance, including standardization, precision, and sensitivity. In certain instances, sensitivity levels of 1 – 10 ng/L is adequate for research and diagnostic applications, however, there is an evolving need for subpicogram/mL quantification to meet the challenges of working with small animals during research, for interrogation of low volume retrospective study samples, and to quantify proteins to evaluate biomarker biology while maintaining a broad dynamic range to accommodate disease progression & therapeutic response. Millipore Sigma's propriety single molecule counting (SMC[™]) immunoassay technology helps to overcome these challenges in the translation of protein biomarkers to the clinic, due to improved assay precision and sensitivity, which has enabled the measurement of endogenous biomarker concentrations at femtogram levels and the precise monitoring over time.





biomarkers in complex biological samples

Figure 2: SMC[™] Assay Workflow This figure illustrates the typical bead based assay work flow. A sandwich immunoassay complex is generated on coated magnetic beads. The complex is disassociated from the bead and the eluate is read on the SMCxPRO[™].



Figure 1: SMCxPRO[™]

Fluorescent based platform, enabling the quantification of low abundance biomarkers in different matrices.





The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada MilliporeSigma and the M logo are trademarks of Merck KGaA, Darmstadt, Germany. Durapore, Millipore Express and Viresolve are registered trademarks of Merck KGaA, Darmstadt, Germany. © 2016 EMD Millipore Corporation, Billerica, MA USA. All rights reserved. Lit. No. PBXXXXXXX XXX

Methods





 Add samples 96-well plate Incubate with capture antibody magnetic beads

 Wash to remove antibodies Incubate with fluorescent-labeled detection antibody

antibody • Transfer to neutralized buffer in 384 plate Beads are not transferred into the read plate

SMC[™] Technology

• Fluorescence captured by optical system and photons counted by APD

Results (II)

Immunoassay	Format	Vendor	LLoQ (pg/mL)	Fold Improvement
Original IL-13 ELISA	ELISA	R&D Systems	9.8	-
Erenna Plate assay	Plate	Singulex	1.2	8×
Erenna MP assay	MP	Singulex	0.07	140×

MOA studies. cohorts).



Figure 3: Single Molecule Counting[™] Technology Single molecule counting is performed inside individual wells. A small 5 μ m integration space is illuminated by a laser. The laser

scans the well, single fluorescently labeled molecules are excited as they pass through the interrogation space. Detected signals with peak intensity above the threshold background of fluorescence are counted as single digital events.

Results (I)



Figure 4: Shifting the Paradigm: Previously undetectable, Quantified

New Biomarker discovery * Monitor small changes in analyte levels * Follow PK studies for longer time courses * Micro-dosing studies for improved tox studies



improving precision.



Figure 7: mHTT is Quantifiable in Mutation Carriers, not Control Participants in 2 Independent Cohorts CSF mHTT levels are significantly elevated in manifest HD compared to premanifest HD carriers.

St. Ledger et al., J Immuno Methods, 2009.

Figure 5: IL-13 as Efficacy Marker for Asthma Accurate measurement of circulating iL-13 levels is needed for PK/PK &

Quantification of baseline levels in clinical study population (n=182, 3)

Using 99th% cut-off, quantified all healthy & mild asthmatic subjects. AII < 1 pg/mL - previously not measurable by other methods

Figure 6: Increase [iL-17A] & [iL-17F] Correlate with RA Ultra-high sensitivity improves S/N response of samples, thereby

Combining the results from the two biomarkers allows for disease identification with improved specificity.

Results (III)

Parameters ass: elements Reagents and reference materia

Target range

Dynamic range (lower and upper quantitation limits)

Sensitivity

Curve fitting Selectivity and specific

Parallelism

Dilution linearity

Precision and accurac (analytical)

Relative accuracy/rec (biological)

Robustness (reagent change control)

Sample handling, coll processing, and stor

Documentation

													Dilution	
		Mean								Mean			Corrected	
		[Progerin]			[Progerin]	Spike			Dilution	[Progerin]			[Progerin]	
n	Sample	pg/mL	SD	%CV	Spike pg/mL	Recovery	n	Sample	Factor	pg/mL	SD	%CV	pg/mL	Linearity
1	Buffer Control	-	-	-			2		5	2,261	60.07	13%	11,305	
2	" rh. Spike	1,048	114.53	11%	1,000	105%	2	Sub#154	25	509	6.81	7%	12,728	113%
2	" endog Spike	847	2.33	0%	1,019	83%	2		125	111	3.79	17%	13,834	109%
2	NHP 1	67	16.52	25%			2		5	5,613	6.04	1%	28,063	
2	" rh. Spike	1,086	68.72	6%	1,000	102%	2	Sub#110	25	985	4.68	2%	24,634	88%
2	" endog Spike	940	81.44	9%	1,019	86%	2		125	200	1.99	5%	25,034	102%
2	NHP 2	-	-	-			2		25	1,075	42.06	0.04	26,865	
2	" rh. Spike	1,254	1.20	0%	1,000	125%	2	Sub#145	50	546	12.31	0.02	27,319	102%
2	" endog Spike	1,058	6.56	1%	1,019	104%	2		100	245	6.97	0.03	24,483	90%
2	NHP 3	31	0.01	0%										
2	" rh. Spike	1,391	108.44	8%	1,000	136%								
2	" endog Spike	1,044	64.55	6%	1,019	99%								
2	Sub#110	550	56.56	10%										
2	" rh. Spike	1,920	2.33	0%	1,000	137%								
2	" endog Spike	1,745	118.29	7%	1,019	117%								
2	Sub#145	559	7.23	1%										
2	" rh. Spike	1,993	93.17	5%	1,000	143%								
2	" endog Spike	1,738	51.78	3%	1,019	116%								

Summary

• Using SMC[™] technology we have developed a sensitive assay for the specific quantification of Progerin

- therapeutics.

Table 1: MS Custom Assays & Sample Testing Group follows Standardized Guidelines

Signat

Established in 2005 "Fit for Purpose" Method Development & Validation for Successful Biomarker Measurement

HUHDORC

Table II. Fit-for-Purpose Elements of Biomarker Assay Development and Method Validation

y	Preanalytical and method development ^o	$\mathbb{E}\mathbf{x}$ ploratory method validation ^b	Advanced method validation ^c			
	Consistent and accessible source do due diligence)	Initial characterization Stablity initiated	Well characterized Inventoried Establish stability Establish change control			
	Estimate in biomarker work plan Define expectation of LLOQ and ULOQ	Acquiring data	Establish from incurred samples			
	Determine preliminary assay range with precision profile over target range	Use 3 validation runs	Use at least 6 runs (in-study validation data can be utilized) Establish LLOQ and ULOQ			
	Define minimum detectable range Define requirements of sensitivity (LOD) and LLOQ	Estimate sensitivity Consider LOD <i>vs.</i> LLOQ	Establish sensitivity			
icity	 Choose appropriate calibration model fitting method and tools Reagent specificity from supplier or literature Assess matrix effects and minimize if possible. Determine minimum required dilution (MRD) Sample and substitute matrix 	Confirm choice of calibration model from 3 validation runs Investigate likely sources of interference, including the therapeutic agent	Use 6 validation runs to confirm calibration model Extensive testing of interference and risk recommendation Assessment of biomarker heterogeneity and isoforms			
	N/A	Use incurred samples, if available	Investigate in targeted population Determine maximum tolerable dilution			
	Determine if applicable, as defined in the biomarker plan (test range)	Use spiked samples	Use spiked samples and dilution VS if applicable			
у	Establish expectations early on in biomarker work plan Consider heterogeneity	Use 3 validation runs	Use of total of at least 6 runs (in study validation data can be utilized)			
overy	Establish expectations early on in biomarker work plan	Use spiked incurred samples at multiple concentrations Addition recovery	Use multiple donors			
and	Determine need Consider availability of biological matrix	NA	Establish tolerability on crucial elements			
ection, rage	Establish feasible conditions	Establish short-term and bench top stability Optimize conditions and effects on assay	Establish freeze/thaw and long-term sample stability			
	Biomarker work plan Draft procedures Assess outsourcing options	Appropriate documentation to support the intended use of the data	Appropriate documentation to support the intended use of the data			

Lee et al., Pharmaceutical Research, Vol. 23, 2006

• The recovery results demonstrate the recombinant standard is appropriate for quantification, whereas the linearity results show accurate measurement at various levels of dilution.

The sensitivity allows the use of mouse models for low volume testing & the specificity enables testing the efficacy of potential